
Product Specifications

VWR Part Number: IR2914
ImmunoReagents, Inc. Part Number: GtxHu-055-U

Product: GOAT ANTI-ALPHA FETOPROTEIN, IG FRACTION

Concentration: >5 mg/ml

Amount: 1.0 mg

Conjugate: N/A

Form: Clear, colorless liquid, 0.2 µm filtered

Purification: Proprietary 2-step fractionation procedure

Purity: >75% based on HPLC

Host: Goat

Immunogen: Purified Human Alpha Fetoprotein (AFP)

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Storage: 2-8 °C

Shelf Life: 1 year from date of receipt. Prepare working dilution prior to use and then discard.

Specificity: Human Alphafetoprotein (AFP)

Cross Reactivity: N/A

Country of Origin: Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.

Applications: For use as a standard or control in most immunoassay formats.

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.