



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

Corporate Headquarters

AMRESKO, LLC

28600 Fountain Parkway
Solon, OH 44139 USA
TEL: 440-349-1199 FAX: 440-349-1182
WWW.AMRESKO-INC.COM
Email: info@amresco-inc.com

CERTIFICATE OF ANALYSIS

A.C.E. Sequencing Buffer, 10X Conc.

Code: K249
Batch Number: 0095C478

Chemical Formula:	Not Applicable	Manufacture Date:	12/11/2015
Molecular Weight:	Not applicable	Reassay Date:	12/31/2019
CAS #:	Not Applicable		
Appearance:		Storage:	Room Temperature
Colorless, haze-free liquid		Grade:	BIOTECHNOLOGY GRADE

TEST	SPECIFICATION	ANALYSIS	DISPOSITION
Conductivity @25C	5.78 - 6.72 mS	6.20 mS	PASS
DNase	NONE / Fail	NONE	PASS
pH @25C	7.60 - 7.80	7.67	PASS
RNase	NONE / Fail	NONE	PASS

Spec Set: K249BIO

Analysis may have been rounded to significant digits in specification limits.
Product meets analytical specifications of the grades listed.

Internal ID #: 546

Signature:

Title: QC Manager

Date Printed: 07/06/2016

Page 1 of 1