



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Giardia lamblia Antibody [6273] (FITC)

CATALOG NUMBER: 35-637

Specifications

SPECIES REACTIVITY:	Protozoan
TESTED APPLICATIONS:	ELISA, IF
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Reactive with Giardia cysts. Non-reactive with Cryptosporidium and Microsporidia cysts.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Giardia lamblia antibody is conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product.
PHYSICAL STATE:	Liquid
BUFFER:	Giardia lamblia antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
STORAGE CONDITIONS:	Giardia lamblia antibody can be stored at -20°C, stable for one year.
CLONALITY:	Monoclonal
ISOTYPE:	Mouse IgG3
CONJUGATE:	FITC

Additional Info

ALTERNATE NAMES:	Anti-Giardia lamblia
-------------------------	----------------------

Background

FOR RESEARCH USE ONLY

December 12, 2016