



Entamoeba histolytica Antibody [3555]

CATALOG NUMBER: 35-570

Specifications

SPECIES REACTIVITY:	Protozoan
TESTED APPLICATIONS:	IHC
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Stains E. histolytica in IHC specimens (formalin, paraffin)w/retrieval. Crossreacts with Toxoplasma gondii.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Entamoeba histolytica antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.
PHYSICAL STATE:	Liquid
BUFFER:	Entamoeba histolytica antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
STORAGE CONDITIONS:	Entamoeba histolytica antibody can be stored at -20°C, stable for one year.
CLONALITY:	Monoclonal
ISOTYPE:	Mouse IgM
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	Anti-Entamoeba histolytica
-------------------------	----------------------------

Background

FOR RESEARCH USE ONLY

December 12, 2016