



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

E. coli Antibody

CATALOG NUMBER: 35-161

Specifications

SPECIES REACTIVITY:

TESTED APPLICATIONS:

USER NOTE: Optimal dilutions for each application to be determined by the researcher.

SPECIFICITY: Specific to all "O" and "K" antigenic serotypes of Escherichia coli

IMMUNOGEN: Mixtures of all antigenic serotypes

HOST SPECIES: Rabbit

Properties

PURIFICATION: Purified IgG fraction of antiserum to a purity of >90%.

PHYSICAL STATE: Liquid

BUFFER: E. coli antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

STORAGE CONDITIONS: E. coli antibody can be stored at -20°C, stable for one year.

CLONALITY: Polyclonal

CONJUGATE: Unconjugated

Additional Info

ALTERNATE NAMES: Anti-E. coli

Background

FOR RESEARCH USE ONLY

December 12, 2016