




Fetuin Antibody

CATALOG NUMBER: 51-864



Immunohistochemistry of human liver tissue stained using Fetuin Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC
APPLICATIONS:	Fetuin antibody can be used in immunohistochemistry starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Recognizes human Alpha 2 HS Glycoprotein.
IMMUNOGEN:	Purified Human Alpha-2 HS glycoprotein.
HOST SPECIES:	Sheep

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2 0.05% sodium azide.
STORAGE CONDITIONS:	Fetuin antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	AHSG, A2HS, Alpha-2-Z-globulin, Fetuin, Fetuin-A, AHS, Alpha-2-HS-glycoprotein, Ba-alpha-2-glycoprotein, FETUA, HSGA
ACCESSION NO.:	P02765
PROTEIN GI NO.:	112910
OFFICIAL SYMBOL:	AHSG
GENE ID:	197

Background

