




ICAM3 Antibody [2F8]

CATALOG NUMBER: 51-866



Immunohistochemistry of human spleen tissue stained using ICAM3 Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC, WB
APPLICATIONS:	ICAM3 antibody can be used in Western Blot, immunohistochemistry starting at 20 ug/mL, and immunoprecipitation.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant protein
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	Store ICAM3 antibody at -20 °C or lower. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG2b,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	ICAM3, CD50, CDW50, ICAM-R, CD50 antigen, ICAM-3
ACCESSION NO.:	P32942
PROTEIN GI NO.:	206729872
OFFICIAL SYMBOL:	ICAM3
GENE ID:	3385

Background

