




SOX9 Antibody [3F11]

CATALOG NUMBER: 51-514



Immunohistochemistry of human colon,
myenteric plexus tissue stained using
SOX9 Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IF, IHC, WB
APPLICATIONS:	SOX9 antibody can be used in Western Blot, and immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	SOX9 (NP_000337, 400 a.a. ~ 510 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	Store SOX9 antibody at -20 °C. Aliquot to avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	SOX9, CMD1, CMPD1, SRY-related HMG-box, gene 9, Transcription factor SOX-9
ACCESSION NO.:	P48436
PROTEIN GI NO.:	1351096
OFFICIAL SYMBOL:	SOX9
GENE ID:	6662

Background

