




LDLRAP1 Antibody [4G4-D5]

CATALOG NUMBER: 51-579



Immunohistochemistry of human small intestine, muscularis propria stained using LDLRAP1 Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IF, IHC, WB
APPLICATIONS:	LDLRAP1 antibody can be used in Western Blot starting at 1:1000 - 1:5000, and immunohistochemistry starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	LDLRAP1 (AAH29770, 1 a.a. ~ 264 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	LDLRAP1 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	LDLRAP1, ARH, ARH1, FHCB1, FHCB2, LDL receptor adaptor protein, ARH2
ACCESSION NO.:	Q5SW96
PROTEIN GI NO.:	116241254
OFFICIAL SYMBOL:	LDLRAP1
GENE ID:	26119

FOR RESEARCH USE ONLY

December 13, 2016