



## C2C12 Cell Slide

CATALOG NUMBER: 17-205

### Specifications

<b>TISSUE CELL TYPE:</b>	C2C12
<b>SPECIES:</b>	Mouse
<b>GENERAL DIAGNOSIS:</b>	C2C12 Cell Slide
<b>AVAILABLE LOTS:</b>	human muscle myoblast
<b>TESTED APPLICATIONS:</b>	ICC
<b>APPLICATIONS:</b>	C2C12 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

### Properties

<b>FIXATION:</b>	4% paraformaldehyde
<b>STORAGE CONDITIONS:</b>	C2C12 cell slide can be stored for up to 6 months at 4°C under desiccate conditions.

### Background

C2C12 (muscle; myoblast)
--------------------------

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

**FOR RESEARCH USE ONLY**

December 12, 2016