



prosci-inc.com



**ProSci Incorporated**  
12170 Flint Place  
Poway, CA 92064



**Toll Free:** +1 (888) 513 9525  
**Local:** +1 (858) 513 2638  
**Fax:** +1 (858) 513 2692

techsupport@prosci-inc.com

## A-20 Cell Slide

CATALOG NUMBER: 17-208

### Specifications

<b>TISSUE CELL TYPE:</b>	A-20
<b>SPECIES:</b>	Mouse
<b>GENERAL DIAGNOSIS:</b>	A-20 Cell Slide
<b>AVAILABLE LOTS:</b>	BALB/cAnN mouse, B lymphocyte, reticulum cell sarcoma
<b>TESTED APPLICATIONS:</b>	ICC
<b>APPLICATIONS:</b>	A-20 cell slide is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunohistochemistry can be directly applied.

### Properties

<b>FIXATION:</b>	4% paraformaldehyde
<b>STORAGE CONDITIONS:</b>	A-20 cell slide can be stored for up to 6 months at 4°C under desiccate conditions.

### Background

A-20 (B lymphocyte; reticulum cell sarcoma)

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

**FOR RESEARCH USE ONLY**

December 12, 2016