



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Human Cell Array II

CATALOG NUMBER: 17-602

Specifications

TISSUE CELL TYPE:	HL-60, Jurkat, K-562, MOLT-4, Raji, THP-1
SPECIES:	Human
GENERAL DIAGNOSIS:	Human Cell Array II
TESTED APPLICATIONS:	ICC
APPLICATIONS:	Human Cell Array II is provided for immunocytochemical analysis of protein and protein-protein interaction. Pretreatment with 0.1% TritonX-100 for 1 minute at room temperature before applying immunocytochemistry is strongly suggested.

Properties

FIXATION:	4% paraformaldehyde
STORAGE CONDITIONS:	Human Cell Array II can be stored for up to 3 months at 4°C under desiccate conditions.

Background

Human cell lines (HL-60, Jurkat, K-562, MOLT-4, Raji and THP-1)

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

FOR RESEARCH USE ONLY

December 12, 2016