



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Mouse Cell Array

CATALOG NUMBER: 17-701

Specifications

| | |
|-----------------------------|--|
| TISSUE CELL TYPE: | A-20, EL4, L1210, NIH/3T3, P815 and RAW 264.7 |
| SPECIES: | Mouse |
| GENERAL DIAGNOSIS: | Mouse Cell Array |
| TESTED APPLICATIONS: | ICC |
| APPLICATIONS: | Mouse Cell Array is provided for immunocytochemical analysis of protein and protein-protein interaction. Procedures for immunocytochemistry can be directly applied. |

Properties

| | |
|----------------------------|--|
| FIXATION: | 4% paraformaldehyde |
| STORAGE CONDITIONS: | Mouse Cell Array can be stored for up to 3 months at 4°C under desiccate conditions. |

Background

Mouse cell lines (A-20, EL4, L1210, NIH/3T3, P815 and RAW 264.7)

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

FOR RESEARCH USE ONLY

December 12, 2016