




## CD45RO Antibody [UCHL1] (FITC)

CATALOG NUMBER: 76-487



Human peripheral blood lymphocytes were stained with FITC UCHL1 with relevant isotype control in Red.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	FACS
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>SPECIFICITY:</b>	The UCHL1 monoclonal antibody specifically reacts with human CD45RO, a 180 kDa isoform of the leukocyte common antigen CD45.
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	The monoclonal antibody was purified utilizing affinity chromatography and unreacted dye was removed from the product.
<b>PHYSICAL STATE:</b>	liquid
<b>BUFFER:</b>	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH7.2.
<b>CONCENTRATION:</b>	5 uL (0.5 ug) / test
<b>STORAGE CONDITIONS:</b>	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	Mouse IgG2a, kappa
<b>CONJUGATE:</b>	FITC

### Additional Info

<b>ALTERNATE NAMES:</b>	LCA, LY5, B220, CD45, L-CA, T200, CD45R, GP180, PTPRC
<b>OFFICIAL SYMBOL:</b>	PTPRC
<b>GENE ID:</b>	5788

### Background

<b>BACKGROUND:</b>	The UCHL1 monoclonal antibody specifically reacts with human CD45RO, a 180 kDa isoform of the leukocyte common antigen CD45. CD45RO is a transmembrane glycoprotein with tyrosine phosphatase activity and is
--------------------	---

expressed by majority of thymocytes, monocytes, granulocytes, and activated memory T lymphocytes. The subsets of peripheral T lymphocytes can be discriminated by using the CD45RO and CD45RA expressing cells.

---

**REFERENCES:**

- 1) Knapp W(1989) Leucocyte typing IV: white cell differentiation antigens. Oxford University Press, 1989.
  - 2) Akbar, A. N., Terry, L., Timms, A., Beverley, P. C., Janossy, G. (1988). Loss of CD45R and gain of UCHL1 reactivity is a feature of primed T cells. The Journal of Immunology, 140(7), 2171-2178.
  - 3) Smith, S. H., Brown, M. H., Rowe, D., Callard, R. E., Beverley, P. C. (1986). Functional subsets of human helper-inducer cells defined by a new monoclonal antibody, UCHL1. Immunology, 58(1), 63.
- 

**FOR RESEARCH USE ONLY**

December 13, 2016