Catalog #:AS-81006

AnaSpec, Inc.

Analytical Data Sheet

Product Name:	5(6)-FAM,SE
Lot No.:	AA26965
Weight:	25mg
Formula:	$C_{25}H_{15}NO_9$
Molecular Weight:	473.4
Appearance:	Orange Powder
Purity:	N/A* *Due to isomerization
Solubility:	1 mmol in 1ml DMF

Remarks: Not for Human Use. Research Purposes Only.

Release By: _____

_____ Date: ____

Quality Control