

Name of the Product : **Lactic Phage Agar**
 Code No. : **M967**
 Section 1 : **Chemical Identification**
 Code No. : M967
 Name of the Product : Lactic Phage Agar
 Produced by : HiMedia Laboratories Pvt. Ltd.
 Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India.
 Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468

Section 2	: Composition			
	Ingredients	Grams/Litre	Ingredients	Grams/Litre
	Casein enzymic hydrolysate	10.00	Yeast extract	5.00
	Beef extract	5.00	Lactose	10.00
	Dipotassium phosphate	5.00	Agar	15.00
	Final pH (at 25°C) 6.8±0.2			

Section 3 : **Hazards Identification**
 Hazard : Not classified as hazardous.

Section 4 : **First - Aid Measures**
 No specific measures necessary.

Section 5 : **Fire Fighting Measures**
 Not combustible.

Section 6 : **Accidental Release Measures**
 No specific measures necessary.

Section 7 : **Handling and Storage**
 Handling - Refer to Section 8
 Storage - Store below 30°C

Section 8 : **Exposure Controls / Personal Protection**
 Wear appropriate Niosh/Msha-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.

Section 9 : **Physical and Chemical Properties**
 Appearance : Homogeneous powder
 Colour : Yellow

Section 10 : **Stability and Reactivity**
 Stability : Product is stable if stored as per the conditions specified under storage of Section No. 7.
 Product loses its potency/performance above 45°C.
 Conditions to avoid : Heat and light.

Section 11 : Toxicological Information

Non toxic.

Section 12 : Ecological Information

Data not available

Section 13 : Disposal Considerations

No special disposal method required except that it be in accordance with current and local authority regulation.

Section 14 : Transport Information

UN No. : Not applicable.

Section 15 : Regulatory Information

Risk Phrases : Not applicable

Safety Phrases: Not applicable

Section 16 : Other Information

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for user's obligation to investigate and determine the suitability of the information for their particular purpose.

