



Item Number AM110
Item Amido Black 10B
CAS Number 1064-48-8
Molecular Formula C₂₂H₁₄N₆O₉S₂Na₂
Molecular Weight 616.52
Absorption Max. (in Water)..... 618 nm

Spectrum Chemicals & Laboratory Products

<i>Corporate Headquarters:</i>	<i>East Coast Facility:</i>
14422 S. San Pedro St.	755 Jersey Ave.
Gardena, CA 90248	New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843	