



Item Number	C1620
Item	alpha-Cyclodextrin
CAS Number	10016-20-3
Molecular Formula	C ₃₆ H ₆₀ O ₃₀
Molecular Weight	972.86
Assay.....	Min. 98%
MAXIMUM LIMITS	
Water.....	10%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843