



Item Number G1024
Item Glycylglycine, Ultrapure
CAS Number 556-50-3
Molecular Formula $C_4H_8N_2O_3$
Molecular Weight 132.12

Assay..... Min. 99.0%
pKa @ 20°C..... 8.0 ± 0.4
pH of a 10% Aqueous Solution..... 5-7
Solubility (10% in Water)..... Clear and colorless

MAXIMUM LIMITS

Loss on Drying..... 0.2%
Residue on Ignition..... 0.2%
Chloride (Cl)..... 0.025%
Sulfate (SO_4)..... 0.03%
Iron (Fe)..... 0.003%
Heavy Metals (as Pb)..... <0.002%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: ***East Coast Facility:***

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843