

**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**1 Identification of the substance/mixture and of the company/undertaking**

- **Product Identifier**
- **Product Name:** TT Broth Base Hajna, Dehydrated
- **Catalog Number:** 249120
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
BD Diagnostic Systems  
7 Loveton Circle  
Sparks, MD 21152  
Telephone.: (410) 771 - 0100 or (800) 638 – 8663
- **Information Department:** Technical Service
- **Emergency telephone number:**  
In case of a chemical emergency, spill, fire, exposure, or accident contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.


**2 Composition/information on ingredients**

- **Chemical characterization:** *Mixture*
- **Description:** *Mixture consisting of the following components.*

**Dangerous Components:**

7647-14-5	sodium chloride	5.43%
-----------	-----------------	-------

**3 Hazards identification**

- **Classification of the substance or mixture**
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**  
 Xi; Irritant  
R36/37/38: Irritating to eyes, respiratory system and skin.
- **Classification system:**  
The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.

**Label elements**

- **Labelling according to EU guidelines:**  
The product has been classified and marked in accordance with regulations on hazardous materials.

**Code letter and hazard designation of product:** Xi Irritant

- **Hazard-determining components of labelling:**  
sodium thiosulfate

(Contd. on page 2)

USA

**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: TT Broth Base Hajna, Dehydrated**

(Contd. of page 1)

**· Risk phrases:**

36/37/38 Irritating to eyes, respiratory system and skin.

**· Safety phrases:**

24/25 Avoid contact with skin and eyes.

22 Do not breathe dust.

36 Wear suitable protective clothing.

7 Keep container tightly closed.

**· NFPA ratings (scale 0-4)**

Health = 1

Fire = 0

Reactivity = 0

**· HMIS ratings (scale 0-4)**

Health = 1

Fire = 0

Reactivity = 0

**4 First aid measures**

- **After inhalation** Supply fresh air; consult doctor in case of complaints.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact**  
Rinse opened eye for 15 minutes under running water. Then consult a doctor.
- **After swallowing** Immediately call a doctor.
- **Information for doctor** Show this product label or this MSDS.

**5 Firefighting measures**

- **Suitable extinguishing agents**  
CO<sub>2</sub>, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Protective equipment:** No special measures required.

**6 Accidental release measures**

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Wipe up with damp sponge or mop.
- **Methods and material for containment and cleaning up:** No special measures required.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.

(Contd. on page 3)

USA

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: TT Broth Base Hajna, Dehydrated**

See Section 13 for disposal information.

(Contd. of page 2)

## 7 Handling and storage

- **Handling**
- **Precautions for safe handling** Ensure good ventilation/exhaustion at the workplace.
- **Information about protection against explosions and fires:**  
No special measures required.
- **Storage**
- **Requirements to be met by storerooms and receptacles:** < 30 °C
- **Information about storage in one common storage facility:**  
Do not store together with oxidizing and acidic materials.  
Store away from metals.
- **Further information about storage conditions:**  
Store in cool, dry conditions in well sealed containers.

## 8 Exposure controls/personal protection

- **Additional information about design of technical systems:**  
No further data; see Section 7.

- **Components with limit values that require monitoring at the workplace:**

### 1317-65-3 calcium carbonate

PEL	15*; 5** mg/m <sup>3</sup> *Total dust **Respirable fraction
REL	10*; 5** mg/m <sup>3</sup> *Total dust **Respirable fraction
TLV	10 mg/m <sup>3</sup> (e)

- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal Protective Equipment**
- **General protective and hygienic measures**  
Keep away from foodstuffs, beverages and feed.  
Immediately remove all soiled and contaminated clothing.  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes and skin.
- **Breathing equipment:**  
In case of brief exposure, use a chemical fume hood or a NIOSH/MSHA-approved respirator.
- **Protection of hands:**



Chemical resistant gloves (i.e. nitrile, or equivalent).

- **Eye protection:** Safety glasses

(Contd. on page 4)

**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: TT Broth Base Hajna, Dehydrated**· **Body protection:** Protective work clothing (lab coat).

(Contd. of page 3)

**9 Physical and chemical properties**· **General Information**· **Appearance:**

<b>Form:</b>	Solid. Powder
<b>Color:</b>	Light green
<b>Odor:</b>	Characteristic

· **pH-value:** 7.6 +/- 0.2· **Change in condition****Melting point/Melting range:** Not determined**Boiling point/Boiling range:** Not determined· **Flash point:** Not applicable· **Auto igniting:** Product is not self igniting.· **Danger of explosion:** Product does not present an explosion hazard.· **Density:** Not determined· **Solubility in / Miscibility with****Water:** Insoluble· **Solvent content:****Solids content:** 100.0 %**10 Stability and reactivity**· **Thermal decomposition / conditions to be avoided:**

No decomposition if used according to specifications.

· **Incompatible materials:** Incompatible with strong acids, oxidizers and lead.· **Hazardous decomposition products:**Carbon monoxide (CO) and carbon dioxide (CO<sub>2</sub>)Sulfur oxides (SO<sub>x</sub>)**11 Toxicological information**· **Acute toxicity:**· **Primary irritant effect:**· **on the skin:** Irritant to skin and mucous membranes.· **on the eye:** Irritating effect.· **Sensitization:** No sensitizing effects known.

(Contd. on page 5)



# Material Safety Data Sheet

acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: TT Broth Base Hajna, Dehydrated**

(Contd. of page 4)

**· Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

## 12 Ecological information

**· Aquatic toxicity:** No further relevant information available.**· Ecotoxicological effects:****· Other information:**

The ecological effects have not been thoroughly investigated, but currently none have been identified.

**· Additional ecological information:****· General notes:** Not known to be hazardous to water.

## 13 Disposal considerations

**· Waste treatment methods****· Recommendation**

Smaller quantities can be disposed of with solid waste.

This product is not considered a RCRA hazardous waste.

Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.

**· Uncleaned packagings:****· Recommendation:** Disposal must be made according to state and federal regulations.**· Recommended cleansing agent:** Water, if necessary with cleansing agents.

## 14 Transport information

**· DOT regulations:****· Hazard class:** -**· Land transport ADR/RID (cross-border)****· ADR/RID class:** -**· Maritime transport IMDG:****· IMDG Class:** -**· Marine pollutant:** No**· Air transport ICAO-TI and IATA-DGR:****· ICAO/IATA Class:** -

USA

(Contd. on page 6)

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: TT Broth Base Hajna, Dehydrated**

(Contd. of page 5)

## 15 Regulatory information

### · SARA Section 355 (extremely hazardous substances)

None of the ingredients is listed.

### · SARA Section 313 (specific toxic chemical listings)

None of the ingredients is listed.

### · TSCA (Toxic Substances Control Act)

7772-98-7 sodium thiosulfate

1317-65-3 calcium carbonate

9000-71-9 casein

7647-14-5 sodium chloride

69-65-8 mannitol

8013-01-2 yeast extract

138-15-8 glutamic acid hydrochloride

50-99-7 dextrose

302-95-4 sodium desoxycholate

9002-18-0 agar

### · California Proposition 65 - Chemicals known to cause cancer

None of the ingredients is listed.

### · California Proposition 65 - Chemicals known to cause reproductive toxicity for females:

None of the ingredients is listed.

### · California Proposition 65 - Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

### · California Proposition 65 - Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

### · Carcinogenic categories

#### · NTP (National Toxicology Program)

None of the ingredients is listed.

#### · TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

### · Product related hazard information:

The product has been classified and marked in accordance with regulations on hazardous materials.

(Contd. on page 7)

USA

**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 02/11/2011

**Product Name: TT Broth Base Hajna, Dehydrated**

(Contd. of page 6)

**· Hazard symbols:**

Xi Irritant

**· Hazard-determining components of labeling:**

sodium thiosulfate

**· Risk phrases:**

36/37/38 Irritating to eyes, respiratory system and skin.

**· Safety phrases:**

24/25 Avoid contact with skin and eyes.

22 Do not breathe dust.

36 Wear suitable protective clothing.

7 Keep container tightly closed.

**16 Other information**

To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.

**· Department issuing MSDS:**

Environmental, Health &amp; Safety

Created by Michael J. Spinazzola

**· Contact:** Technical Service Representative

— USA —