

NBP1-91284**IFN gamma Receptor beta Protein****Unit Size**

0.11 mg

Concentration

Lyoph

Storage

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

Applications

ELISA, FLOW

Species

Human

Purity

>95% pure by SDS-PAGE

Specificity/Sensitivity

33 kDa dimer composed of identical 16.5 kDa monomer Biological activity has been assessed in a bioassay using inhibition of the cytopathic encephalomyocarditis virus effect on human lung carcinoma cell line as test parameter.

Recommended Dilutions

ELISA

Flow Cytometry

Buffer

PBS containing 10 mg/ml bovine serum albumin. Reconstitute with 1.1 ml distilled water

Uses

This protein is useful for ELISA and Flow Cytometry

ImagesThis product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: www.novusbio.com/NBP1-91284

Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

