

**NBP1-76926****Nucleoporin 107 Antibody****Unit Size**

0.1 mg

**Concentration**

1 mg/ml

**Storage**

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

**Applications**

WB, ELISA, ICC/IF

**Species**

Human, Mouse, Rat

**Host**

Rabbit

**Clonality**

Polyclonal

**Purity**

Affinity purified

**Preservative**

0.02% Sodium Azide

**Immunogen**

Antibody was raised against a 20 amino acid peptide from near the center of human NUP107.

**Recommended Dilutions**

ELISA, Immunocytochemistry/Immunofluorescence, Western Blot

**Buffer**


PBS containing 0.02% sodium azide.

**Uses**

ELISA, Western Blot, Immunocytochemistry

**Images**

Western Blot: Nucleoporin 107 Antibody [NBP1-76926] - A549 cell lysate with NUP107 antibody at (A) 1 and (B) 2 ug/mL.



Immunocytochemistry/Immunofluorescence:  
 Nucleoporin 107 Antibody [NBP1-76926] - A549 cells with Nucleoporin 107 antibody at 10 ug/ml.

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

