

NBP1-27972**clAP2 Antibody****Unit Size**

0.1 ml

Concentration

1.5 mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Applications

WB, ICC/IF, IHC, IHC-P

Species

Human

Molecular Weight

64kDa

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified

Preservative

0.1% Sodium Azide

Immunogen

Synthetic peptide made to an internal portion of human clAP2 (within residues 100-200). [Swiss-Prot# Q13489]

Recommended Dilutions

Immunocytochemistry/Immunofluorescence 1:100, Immunohistochemistry 1:800-1:3200, Immunohistochemistry-Paraffin 1:800-1:3200, Western Blot 1.0 ug/ml

Buffer

PBS, 30% glycerol

Uses


This clAP2 antibody is useful for Immunocytochemistry/Immunofluorescence, Immunohistochemistry paraffin embedded sections and Western blot analysis, where a band can be seen at ~64 kDa.

References (more available at www.novusbio.com/NBP1-27972)


Guicciardi ME, Mott JL, Bronk SF et al. Cellular Inhibitor of Apoptosis 1 (clAP-1) Degradation by Caspase 8 During TNF-related Apoptosis-inducing Ligand (TRAIL)-induced Apoptosis. Exp Cell Res. 2010 Oct 13. [PMID: 20951133]

Images (more available at www.novusbio.com/NBP1-27972)

Western Blot: clAP2 Antibody [NBP1-27972] - HeLa whole cell lysate.



Immunohistochemistry: clAP2 Antibody [NBP1-27972] - Analysis of human adenocarcinoma cells and adjacent glands using anti-clAP2.



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA
 8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com
Novus Canada
 461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com
Novus Europe
 12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com
