

NBP1-06652

## NKG2D Antibody (AP-MAB0840)

### Unit Size

0.5 mg

### Concentration

1.0 mg/ml

### Storage

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

### Applications

FLOW

### Species

Mouse

### Host

Armenian Hamster

### Clonality

Monoclonal

### Isotype

IgG2

### Purity

Protein G purified

### Preservative

No Preservative

### Specificity/Sensitivity

This was produced from a hybridoma (mouse myeloma fused with spleen cells from a Hamster immunized with Soluble mouse NKG2D).

### Immunogen

Soluble mouse NKG2D.

### Recommended Dilutions

Flow Cytometry ,

### Buffer

0.5mg antibody in 0.5 mL PBS with no preservative

### Uses

This antibody is useful for Flow Cytometry as well as Blocking NK cell killing of targets expressing NKG2D ligands. Extremely low level of LPS (&lt; 0.002EU/ug IgG)

### Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/NBP1-06652](http://www.novusbio.com/NBP1-06652)

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

#### Novus USA

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

#### Novus Canada

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

#### Novus Europe

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

