

NB600-434

## SMC3 Antibody [Agarose Immobilized]

### Unit Size

0.1 mg

### Concentration

Please see the vial label for concentration.

### Storage

Store at 4C. Do not freeze.

### Applications

IP

### Species

Human, Xenopus

### Host

Rabbit

### Clonality

Polyclonal

### Conjugate

Agarose Immobilized

### Purity

Immunogen affinity purified

### Preservative

0.09% Sodium Azide

### Immunogen

Rabbit anti-human SMC3 affinity purified Antibodies were coupled to agarose beads using a cyanogen bromide method.

### Recommended Dilutions

Immunoprecipitation

### Buffer

Phosphate Buffered Saline (PBS)

### Uses

This antibody is useful for Immunoprecipitation 15 -25 mcl of gel slurry per 1 -1 mg of protein lysate or extract.

### Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/NB600-434](http://www.novusbio.com/NB600-434)

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

#### Novus USA

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

#### Novus Canada

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

#### Novus Europe

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

