

NB400-163B**ABCA7 Antibody (7A1-144) [Biotin]****Unit Size**

0.1 ml

Concentration

1.3 mg/ml

Storage

Store at 4C in the dark.

Applications

WB, FLOW

Species

Mouse

Host

Rat

Clonality

Monoclonal

Isotype

IgG2a Kappa

Conjugate

Biotin

Purity

Protein G purified

Preservative

No Preservative

Specificity/Sensitivity

This is specific for ABCA7.

Immunogen

Mouse ABCA7 stably-transfected HeLa cells.

Recommended Dilutions

Flow Cytometry 1:500 - 1:4000 , Western Blot 1:2000

Buffer

PBS

Uses

By Western Blot a 240 kDa band is seen. NB 400-163 is not reactive in ELISA.

Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: www.novusbio.com/NB400-163B

Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. For more information on our guarantee, please visit www.novusbio.com/guarantee.

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

