

**NB110-75118****OR11H4 Antibody****Unit Size**

0.1 ml

**Concentration**

lyoph

**Storage**

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

**Applications**

WB, IHC, IHC-Fr

**Species**

Human

**Host**

Rabbit

**Clonality**

Polyclonal

**Purity**

Whole antisera

**Preservative**

No Preservative

**Specificity/Sensitivity**

Human olfactory receptor 11H4.

**Immunogen**

A synthetic peptide as a part of human olfactory receptor 11H4 conjugated to immunogenic carrier protein was used as the immunogen. The antigen is homologous with the corresponding sequence derived from human olfactory receptors OR11G2 and OR14-34.

**Recommended Dilutions**

Western Blot 1:400 to 1:3000, Immunohistochemistry-Frozen, Immunohistochemistry

**Buffer**

Lyophilized Reconstitute in 0.1 ml of sterile water. Centrifuge to remove any insoluble material. When reconstituting, glycerol (1:1) may be added for an additional stability.

**Uses**

Immunohistochemistry-Paraffin/Frozen and Western Blot

**Images**This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/NB110-75118](http://www.novusbio.com/NB110-75118)

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

