

NB100-65058**Shigella Group Antigen Antibody****Unit Size**

0.5 ml

Concentration

5.0 mg/ml

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Applications

WB, ELISA, ICC/IF, IHC, IHC-Fr

Species

Bacteria

Host

Rabbit

Clonality

Polyclonal

Purity

IgG purified

Preservative

0.09% Sodium Azide

Specificity/Sensitivity

NB100-65058 is polyvalent for Shigella antigens. *Shigella sp.* are non spore forming, rod shaped bacteria, the causative agent of Shigellosis in man. The genus is closely related to *Escherichia coli* and *Salmonella*. The antibody preparation has not been cross absorbed and is likely to react with closely related Enterobacteriaceae

ImmunogenMixture of *S. boydii*, *S. flexneri* and *S. dysenteriae*.**Recommended Dilutions**

Immunohistochemistry-Frozen, ELISA, Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry

Buffer

PBS

Uses

Immunohistochemistry-Frozen, ELISA, Immunofluorescence and Western Blot

Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: www.novusbio.com/NB100-65058

Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

