

NB100-56493**CRM1 Antibody****Unit Size**

0.1 mg

Concentration

0.5 mg/ml

Storage

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

Applications

WB, IHC, IHC-P

Species

Human, Bovine, Canine, Primate

Molecular Weight

123kDa

Host

Rabbit

Clonality

Polyclonal

Purity

Protein G purified

Preservative

0.05% Sodium Azide

Specificity/Sensitivity

The amino acid sequence used as immunogen has 100% identity in human, chimp (isoforms 1-7), bovine, chicken and dog, and 94% identity in mouse and rat.

Immunogen

Amino acids 390-408 (ASPLLSGSQHFDVPPRRQL) of human Exportin-1 protein were used as the immunogen (NP_003391).

Recommended Dilutions

Immunohistochemistry 10 ug/ml

Immunohistochemistry-Paraffin 10 ug/ml


Western Blot 1-3 ug/ml

Buffer


PBS containing 0.2% gelatin

Images

Western Blot: Exportin-1 Antibody [NB100-56493] - Western blot analysis of Exportin-1 in Jurkat cell lysate in the 1) absence and 2) presence of immunizing peptide using NB100-56493 at 1.5 ug/ml.



Immunohistochemistry: Exportin-1 Antibody [NB100-56493] - IHC analysis of human spleen using NB100-56493 at 10 ug/ml.



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

