

NB100-359H**ATRIP Antibody [HRP]****Unit Size**

0.1 ml

Concentration

1.0 mg/ml

Storage

Store at 4C in the dark.

Applications

WB, IP

Species

Human

Host

Rabbit

Clonality

Polyclonal

Conjugate

HRP

Purity

Protein G purified

Preservative

No Preservative

Immunogen

A synthetic peptide derived from a sequence between amino acids 500-550 of the human ATRIP protein (aka TREX1) [Gene ID 11277].

Recommended Dilutions

immunoprecipitation 3-8 ug/ml

Western Blot 1:5000-1:15000

Buffer

PBS

Uses

This antibody may be used for Western blot and immunoprecipitation. Immunohisto/cytochemistry have not been tested using this antibody.

Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: www.novusbio.com/NB100-359H

Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

