

## Datasheet

### APLP2 purified MaxPab rabbit polyclonal antibody (D01P)

**Catalog Number:** H00000334-D01P

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against a full-length human APLP2 protein.

**Immunogen:** APLP2 (AAH04371.1, 1 a.a. ~ 522 a.a) full-length human protein.

**Sequence:**

MAATGTAAAAATGRLLLLLVGLTAPALALAGYIEALAA  
NAGTGFAVAEPQIAMFCGKLNMHVNIQTGKWEPTPT  
GTKSCFETKEEVLQYCQEMYPELQITNVMEANQRVSI  
DNWCRRDKKQCKSRFVTPFKCLVPPTPLPTNDVDVYF  
ETSADDNEHARFQKAKEQLEIRHRNRMDRVKKEWEE  
AELQAKNLPKAERQTLIQHFQAMVKALEKEAASEKQQ  
LVETHLARVEAMLNDRRRMALENYLAALQSDPPRPHR  
ILQALRRYVRAENKDRDLHTIRHYQHVLAVDPEKAAQMK  
SQVMTHLHVIEERRNQTLSELLKVPYVAQEIQEEIDELL  
QEQRADMDQFTASISETPVDVRVSSEEEIIPPFHPFH  
PPFALPENEGSGVGEQDGLIGAEKVINSKNKVDEN  
MVIDETLDVKEMIFNAERVGGLEERESVGPLREDFSL  
SSSALIGLLVIAVAIATVIVISLVMLRKRQYGTISHGIVEV  
DPMLTPEERHLNKMQNHGYENPTYKYLEQMCI

**Host:** Rabbit

**Reactivity:** Human

**Applications:** WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 334

**Gene Symbol:** APLP2

**Gene Alias:** APPH, APPL2, CDEBP

**References:**

1. Abnormal Platelet Amyloid- $\beta$  Protein Precursor (A $\beta$ PP) Metabolism in Alzheimer's Disease: Identification and Characterization of a New A $\beta$ PP Isoform as Potential Biomarker. Jelic V, Hagman G, Yamamoto NG, Teranishi Y, Nishimura T, Winblad B, Pavlov PF J Alzheimers Dis. 2013 Feb 6.