

Datasheet

ARAF purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00000369-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human ARAF protein.

Immunogen: ARAF (NP_001645.1, 1 a.a. ~ 606 a.a) full-length human protein.

Sequence:

MEPPRGPPANGAEPsRAVGTVKVYLPNKQRTVVTVR
DGMSVYDSLKDALKVRGLNQDCCVVYRLIKGRKTVTA
WDTAIAPLDGEELIVEVLEDVPLTMHNFVRKTFSLAF
CDFCLKFLFHGFRCQTCGYKFHQHCSSKVPTVCVDM
STNRQQFYHSVQDLGGSRQHEAPSNRPLNELLTPQ
GPSRPTQHCDPEHFPPAPANAPLQRIRSTSTPNVHM
VSTTAPMDSNLIQLTGQSFSTDAAGSRGSDGTGRGS
PSPASVSSGRKSPHSKSPAEQRERKSLADDDKKVKNL
GYRDSGYYWEVPPSEVQLLKRIGHTGSFGTVFRGRWH
GDVAVKVLKVSQPTAEQAQAFKNEMQVLRKTRHVNIL
LFMGFMTRPGFAITQWCEGSSLYHHLHVADTRFDMV
QLIDVARQTAQGM DY LHAKNIIHRDLKSNNIFLHEGLTV
KIGDFGLATVKTRWGAQPLEQPSGSLVWMAAEVIRM
QDPNPYSFQSDVYAYGVVLYELMTGSLPYSHIGCRDQ
IIFMVGRGYLSPDLSKISSNCPKAMRRLLSDCLKFQRE
ERPLFPQILATIELLQRLPKIERSASEPSLHRTQADEL
PACLLSAARLVP

Host: Rabbit

Reactivity: Human

Applications: IF, PLA-Ce, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 369

Gene Symbol: ARAF

Gene Alias: A-RAF, ARAF1, PKS2, RAFA1