

Datasheet

RHOB purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00000388-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human RHOB protein.

Immunogen: RHOB (NP_004031.1, 1 a.a. ~ 196 a.a) full-length human protein.

Sequence:

MAAIRKKLVVVGDGACGKTCLLIVFSKDEFPEVYVPTV
FENYVADIEVDGKQVELALWDTAGQEDYDRLRPLSYP
DTDVILMCFSDSPDSLENIPEKWVPEVKHFCPNVPIL
VANKKDLRSDEHVRTELARMKQEPVRTDDGRAMAVRI
QAYDYLECSAKTKEGVREVFETATRAALQKRYGSQN
GCINCKVL

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 388

Gene Symbol: RHOB

Gene Alias: ARH6, ARHB, MST081, MSTP081, RHOH6