

Datasheet

ARR3 monoclonal antibody (M07), clone 2D7

Gene Alias: ARR3

Catalog Number: H00000407-M07

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant ARR3.

Clone Name: 2D7

Immunogen: ARR3 (AAH12096, 1 a.a. ~ 359 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MSKVFKKTSSNGKLSIYLGKRDFVDHVDTVEPIDGVVL
VDPEYFKCRKLFVMLTCAFRYGRDDLEVIGLTFRKDL
VQTLQVVPAAESSSPQGPLTVLQERLLHKLGDNAYPFT
LQMVNTNLPCSVTLQPGPEDAGKPCGIDFEVKSFCAEN
PEETVSKRDYVRLVVRKVQFAPPEAGPGPSAQTIRRF
LLSAQPLQLQAWMDREVHYHGEPISVNVSINNCTNKVI
KKIKISVDQITDVVLYSLDKYTKTVFIQEFTETVAANSS
SQSFAVTPILAASCQKRGLALDGKLKHEDTNLASSTIIR
PGMDKELLGILVSYKVRVNLMSVSCGILGDLTASDVGV
ELPLVLIHPKPSYEAAR

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 407

Gene Symbol: ARR3