

Datasheet

PASD1 monoclonal antibody (M07), clone 1F5

Catalog Number: H00139135-M07

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant PASD1.

Clone Name: 1F5

Immunogen: PASD1 (NP_775764, 1 a.a. ~ 100 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MKMRGEKRRDKVNP KSSQRKLNWIPSFPTYDYFNQV
TLQLLDGFMITLSTDGVII CVAENISSLLGHLPAEIVGKK
LLSLLPDEEKDEVYQKIILKFPLL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 139135

Gene Symbol: PASD1

Gene Alias: -