

## Datasheet

### UPRT purified MaxPab mouse polyclonal antibody (B01P)

**Catalog Number:** H00139596-B01P

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length human UPRT protein.

**Immunogen:** UPRT (NP\_659489.1, 1 a.a. ~ 309 a.a) full-length human protein.

**Sequence:**

MATELQCPDSMPCHNQVNSASTPSPEQLRPGDLILD  
HAGGNRASRAKVILLTGYAHSSSLPAELDSGACGGSSL  
NSEGNSGSGDSSSYDAPAGNSFLEDCELSRQIGAQLK  
LLPMNDQIRELQTIIRDKTASRGDFMFSADRLIRLVVEE  
GLNQLPYKECMVTTPTGYKYEGVKFEKGNGVVSIMRS  
GEAMEQGLRDCCRSIRIGKILIQSDEETQRAKVYYAKF  
PPDIYRRKVLLMYPILSTGNTVIEAVKVLIEHGVQPSVIL  
LSLFSTPHGAKSIIQEFPEITILTTEVHPVAPTHFGQKYF  
GTD

**Host:** Mouse

**Reactivity:** Human

**Applications:** Det Ab, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 139596

**Gene Symbol:** UPRT

**Gene Alias:** DKFZp781E1243, FUR1, MGC23937, UPP