

Datasheet

C20orf132 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00140699-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human C20orf132 protein.

Immunogen: C20orf132 (AAH57767.1, 1 a.a. ~ 608 a.a) full-length human protein.

Sequence:

MPLSRRGRGSRLAGSRAATSNLPRATALSPKSAGRG
ALSISHAQLLTRRRHSSEQVPPGSEPRADFRSGKWFO
EPATGDARDSRQALRARMSSKHRICSQEEVVIPCAVD
SDSESVDLELSNLEIHKKGSSSIELTDLDIPDIPGLHCEP
LSHSPRHILTQQDPLSEAIVEKLIQSIQKVFNGELKGELE
KLKFLGDLSSLSQALPYDETAKSFIHSHIADIVHTLNVLV
QEERPHSLSSSMRQEVFVTIADLSYQDVHLLGSEDR
AELFSLTIKSIITLPSVRTLTQIQEIMPNGTCNTECLYRQ
TFQAFSEMLQSLVVKDPHLENLDTIHKHLPWLQSVKD
HERERATASMAQVLKCLSKHLNLKLPLRFQRLGHLVA
LMALLCGDPQEKVAEEAAEGIHSLHITLRLKYITHDKK
DQQNLKRALTKREFLELHSSAAKCFYNCPFRIAQVFE
GFLDSNELCQFIMTTFTLTKLKHPCIQRSAGELLLTLA
KNTESQFEKVPEIMGVICAQLPIISQPRVRQIINTVSLF
ISRPKYTDIVLSFLLCHPVPYNRHLAEVWRMLSVELPS
TTWILWRLRLKLQKCHNEPAQEKMAYVAVAVSP

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 140699

Gene Symbol: C20orf132

Gene Alias: C20orf131, DKFZp434N0426, FLJ36113, dJ621N11.3, dJ621N11.4