

Datasheet

Cyp2b1 monoclonal antibody

Catalog Number: MAB3815

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against Cyp2b1.

Immunogen: Rat Cyp2b1.

Host: Mouse

Reactivity: Rat

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Identification of 2B1 by immunoassay or Western blots, without cross-reactivity with 2B2 or 2A1.

Form: Liquid

Isotype: IgG

Recommend Usage: Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In purified antibody

Storage Instruction: Store at -20°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 24300

Gene Symbol: Cyp2b1

Gene Alias: Cyp2b-1, Cypbe, MGC188065