

Datasheet

LOC143241 monoclonal antibody (M01), clone 1B11

Catalog Number: H00143241-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant LOC143241.

Clone Name: 1B11

Immunogen: LOC143241 (AAH19250, 1 a.a. ~ 177 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MESIYLQKHLGACLTQGLAEVARVRPVDPIEYLALWIY
KYKENVTMEQLRQKEMAKLERERELALMEQEMMERL
KAEELLQQQQLALQLELEMQEKERQRIQELQRAQEQ
LGKEMRMNMENLVRNEDILHSEEATLDSGKTLAEISDR
YGAPNLSRVEELDEPMFSDIALNIDQDL

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 143241

Gene Symbol: DYDC1

Gene Alias: DPY30D1, FLJ43920, bA36D19.5