

Datasheet

MIPOL1 polyclonal antibody (A01)

Catalog Number: H00145282-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant MIPOL1.

Immunogen: MIPOL1 (NP_620059, 321 a.a. ~ 423 a.a) partial recombinant protein with GST tag.

Sequence:

ARETAVQQYKKLEEEIQTLRVYYSLHKSLSQEENLKDQ
FNYTLSTYEEALKNRENIVSITQQQNEELATQLQQALT
ERANMELQLQHAREASQVANKEVQKLE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 145282

Gene Symbol: MIPOL1

Gene Alias: DKFZp313M2036, MGC34010