

Datasheet

ZNF688 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00146542-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF688 protein.

Immunogen: ZNF688 (NP_660314, 1 a.a. ~ 276 a.a)
full-length human protein.

Sequence:

MAPPPAPLLAPRPGETRPGCRKPGTVSFADVAVYFSP
EEWGCLRPAQRALYRDVMQETYGHLGALGFPGPKPA
LISWMEQESEAWSAAQDPEKGERLGGARRGDVPNR
KEEEPEEVPRAKGPRKAPVKESPEVLVERNPDPAISV
APARAQPPKNAAWDPTTGAQPPAPIPSMDAQAGQRR
HVCTDCGRRFTYPSLLVSHRRMHSGERPFPCECGM
RFKRKFAVEAHQWIHRSCSGGRRGRRPGIRAVPRAP
VRGDRDPPVLFRHYDPDIFEECG

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 146542

Gene Symbol: ZNF688

Gene Alias: -