

Datasheet

WDR16 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00146845-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human WDR16 protein.

Immunogen: WDR16 (AAH25392.1, 1 a.a. ~ 620 a.a) full-length human protein.

Sequence:

MDNKISPEAQVAELELDAVIGFNGHVPTGLKCHPDQE
HMIYPLGCTVLIQAINTEQNFLQGHGNNVSCLAISRS
GEYIASGQVTFMGFKVDIILWDYKNRELLARLSLHKGKI
EALAFSPNDLYLVSLGGPDDGSVVVWSIAKRDAICGSP
AAGLNVGNATNVIFSRCRDEMFM TAGNGTIRVWELDL
PNRKIWPTESQTGQLKRIVMSIGVDDDDSFYLGTTTG
DILKMNPRTKLLTDVGPADKDFSLGVS AIRCLKMGGLL
VGSGAGLLVFCKSPGYKPIKKIQLQGGITSITLRGEGH
QFLVGTEESHIYRVSFTDFKETLIATCHFDAVKDIVFPF
GTAELFATCAKKDIRVWHTSSNRELLRITVPNMTCHGI
DFMRDGKSIISAWNDGKIRAFAPETGRLMYVINNAHRI
GVTAIATTS DCKRVISGGGEGEVRVWQIGCQTQKLEE
ALKEHKSSVSCIRVKRNNEECVTASTDGTCTIWDLVRL
RRNQMILANTLFQCVCYHPEEFQITSGTDRKIAYWEV
FDGTVIRELEGSLSGSINGMDITQEGVHFVTGGNDHLV
KVWDYNEGEVTHVGVGHSGNITRIRISPGNQYIVSVSA
DGAILRWKYPYTS

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 146845

Gene Symbol: WDR16

Gene Alias: FLJ37528, WDRPUH

Gene Summary: WD repeat-containing proteins, such as WDR16, play crucial roles in a wide range of physiologic functions, including signal transduction, RNA processing, remodeling the cytoskeleton, regulation of vesicular traffic, and cell division (Silva et al., 2005 [PubMed 15967112]).[supplied by OMIM]