

Datasheet

C18orf23 monoclonal antibody (M01), clone 3A1

Catalog Number: H00147341-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant C18orf23.

Clone Name: 3A1

Immunogen: C18orf23 (NP_689683, 122 a.a. ~ 220
a.a) partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

VSHPRRSQERVSVHPhRLHPSFDGQLQTPQPRYLA
EGTDWDLSDAGLSPAQFQVRPIQHYQHLYLATPRM
HHFPRNSSSTQMVVHEIRNYPYPQLHF

Host: Mouse

Reactivity: Human,Rat

Applications: ELISA, IF, WB-Re, WB-Ti
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 147341

Gene Symbol: C18orf23

Gene Alias: FLJ34218