

Datasheet

HIPK4 MaxPab mouse polyclonal antibody (B01)

Entrez GeneID: 147746**Gene Symbol:** HIPK4**Gene Alias:** FLJ32818**Catalog Number:** H00147746-B01**Regulation Status:** For research use only (RUO)**Product Description:** Mouse polyclonal antibody raised against a full-length human HIPK4 protein.**Immunogen:** HIPK4 (NP_653286.2, 1 a.a. ~ 616 a.a) full-length human protein.**Sequence:**

MSTIQSETDCYDIEVLGKGTGGEVAKGWRRSTGEMV
AIKILKNDAYRNRIIKNELKLLHCMRGLDPEEAHVIRFLE
FFHDALKFYLVFELLEQNLFEFQKNNFAPLPARHIRT
TLQVLTALARLKLAIHADLKPENIMLVDQTRCPFRVK
VIDFGSASIFSEVRYVKEPYIQSRFYRAPEILLGLPFCEK
VDVWSLGCVMAELHLGWPLYPGNNEYDQVRYICETQ
GLPKPHLLHAACKAHFFKRNPHPDAAANPWQLKSSA
DYLAETKVRPLERRKYMLKSLDQIETVNGGSVASRLTF
PDREALAEHADLKSMVELIKRMLTWESHERISPSAALR
HPFVSMQQLRSAHETTHYYQLSLRSYRLSLQVEGKPP
TPVVAEDGTPYYCLAEKEAAGMGSVAGSSPFFREE
KAPGMQRAIDQLDDLSLQEAGHGLWGETCTNAVSDM
MVPLKAAITGHHVPDSGPEPILAFYSSRLAGRHKARKP
PAGSKSDSNFSNLIRLSQVSPEDDRPCRGSSWEEGE
HLGASAEPLAILQRDEDGPNIDNMTMEAERPDPFLFD
PSSCPGEWLSEPDCTLESVRGPRAQGLPPRRSHQH
PPRGATSFLQHVTGHH

Host: Mouse**Reactivity:** Human**Applications:** IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Storage Buffer:** No additive**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.