

Datasheet

ZNF599 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: ZNF599

Gene Alias: -

Catalog Number: H00148103-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF599 protein.

Immunogen: ZNF599 (NP_001007249.1, 1 a.a. ~ 588 a.a) full-length human protein.

Sequence:

MAAPALALVSFEDVVVFTFTGEEWGHLDLAQRTLYQEV
MLETCRLLVSLGHPVPKPELIYLLHGGQELWTVKRGLS
QSTCAGEKAKPKITEPTASQLAFSEESSFQELLAQRSS
RDSRLGQARDEEKLKIQEGNLRPGTNPHKEICPEKLS
YKHDDLEPDDSLGLRVLQERVTPQDALHECDSQGGPG
KDPMTDARNNPYTCTECGKGFSKKWALVRHQIHAG
VKPYECNECGKACRYMADVIRHMLHTGEKPYKCIEC
GKAFKRRFHLTEHQRIHTGDKPYECKECGKAFTHRSS
FIQHNMTHTREKPFLECKECGKAFYYSSSFAQHMRIHT
GKKLYECGECGKAFTHRSTFIQHNVTHTGEKPFLECKE
CGKTFCLNSSFTQHMRIHTGEKPYECGECGKAFTHRS
TFIRHKRTHHTGEKPFLECKECGKAFCDSLSLIQHMRIHT
GEKPYECSECGKAFTHHSVFIRHNRTHSGQKPLECKE
CAKAFYYSSSFTRHMRIHTGEKPYVCRECGKAFTQPA
NFVRHNRHTGEKPFLECKEKEKAFCDNFALTQHMRTH
TGEKPFECNECGKTFSSHSSSFTHHRKIHTRV

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 148103