

Datasheet

UBE2U polyclonal antibody (A01)

Catalog Number: H00148581-A01

Regulatory Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant UBE2U.

Immunogen: UBE2U (NP_689702.1, 9 a.a. ~ 102 a.a) partial recombinant protein with GST tag.

Sequence:

LHRDFCDLKENNYKGITAKPVSEDMMEWEVEIEGLQN
SVWQGLVFQLTIHFTSEYNYAPPVVKFITIPFHPNVDPH
TGQPCIDFLDNPEKWNTN

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 148581

Gene Symbol: UBE2U

Gene Alias: MGC35130, RP4-636O23.1