

Datasheet

CCT8L2 monoclonal antibody (M09), clone 1G8

Catalog Number: H00150160-M09

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant CCT8L2.

Clone Name: 1G8

Immunogen: CCT8L2 (AAH33797, 1 a.a. ~ 557 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MDSTVPSALELPQRLALNPRESPRSPEEEEPHLLSSLA
AVQTLASVIRPCYGPGRQKFLVTMKGETVCTGCATAI
LRALELEHPAAWLLREAGQTQAENSGDGTAFVLLTE
ALLEQAEQLLKAGLPRPQLREAYATATAEVLATLPSLAI
QSLGPLEDPSWALHSMNTHLSPMDHLTKLVAHAC
WAIKELDGSFKPERVGVCALPGGTLEDSCLLPGLAISG
KLGGQMATVLSGARVALFACPFPAHPNAPATARLSS
PADLAQFSKGSQDQLEKQVQGQLAAAGINVAVVLGEVD
EETLTADKYGIVVIQARSWMEIYLVSEVLDTPLLPRLLP
PQRPGKCQRVYRQELGDGLAVVFEWECTGTPALTVV
LRGATTQGLRSAEQAVYHGIDAYFQLCQDPRLIPGAG
ATEMALAKMLSDKGSRLGPSGPAFLAFWALKYLPK
TLAENAGLAVSDVMAEMSGVHQGNLLMGVGTGEIIN
VAQEGVWDTLIVKAQGFRAVAEVLQLVTVDIIVVAKK
SPTHQEIWNPDSKKTCKHPPPVETKILGLNN

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 150160

Gene Symbol: CCT8L2

Gene Alias: CESK1, MGC118839, MGC118840