

Datasheet

RP5-821D11.2 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00150365-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human RP5-821D11.2 protein.

Immunogen: RP5-821D11.2 (AAH32248, 1 a.a. ~ 238 a.a) full-length human protein.

Sequence:

MLALTLAKADSPRTALLCSAWLLTASFSAQQHKGSLQ
VHQTLSEMDQVLKALSFPKKKAALLSAAILCFLRTALR
QSFSSALVALVPSGAQPLPATKDTVLAPLRMSQVRS
VIGLQNLLVQKDPLLSQACVGCLEALLDYLDARSPDIAL
HVASQPWNRFLFTLLDAGENSFLRPEILRLMTLLQSM
GHLADHSMAQTLQASLEGLPPSTSSGQPPLQDMLCL
GGVAVSLSHIRN

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 150365

Gene Symbol: MEI1

Gene Alias: MGC40042