

Datasheet

ARL6IP6 purified MaxPab mouse polyclonal antibody (B02P)

Catalog Number: H00151188-B02P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ARL6IP6 protein.

Immunogen: ARL6IP6 (NP_689735.1, 1 a.a. ~ 226 a.a) full-length human protein.

Sequence:

MSFAESGWRSALRRRGPGTPGPVARPSYSSFTQGDS
WGEDEVDEEEGCDQVARDLRAEFSAGAWSEPRKRS
VLPPDGNNGSPVLPDKRNGIFPAAAGSRAQPRRWPVQ
VLSILCSLLFAILLAFLLAIAYLIVKELHAENLKNEDDVT
GLLGFWTLLIISLTAGFSCCSFSWTVTYFDSFEPGMFP
PTPLSPARFKKLTGHSFHMGYSMAILNGIVAALTVAWC
LM

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 151188

Gene Symbol: ARL6IP6

Gene Alias: MGC33864, PFAAP1