

Datasheet

ALS2CR11 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00151254-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ALS2CR11 protein.

Immunogen: ALS2CR11 (AAH30659, 1 a.a. ~ 623 a.a) full-length human protein.

Sequence:

MEPPQETNRPFSTLDNRSGQVQVLSATPLLQRNPYSS
PDIMHIKGSEASSVPYALNQGTTALPKNKNQEGTGHR
LLNMLRKTLESDSEELEITQETPNLVPFGDVVGCLGI
HIKNCRHFMPKISLQHYANLFIKISINKAVKCTKMCSLLS
KNDEKNTVIKFDEVKYFSVQVPRRYDDKRNNILLELIQY
DNREKRAFLGSGVQIHLVEVIQKGCFFIEEVQVLHGNIFV
CRLEVEFMFSYGNFGYGFHQLKPLQKITEPSMFMNL
APPPERTDPVTKVITPQTVEYPAFLSPDLNVTGTPAV
QSSNQPSVVRLEKLQQPRERLEKMKKEYRNLNTWI
DKANYLESILMPKLEHKDSEETNIDEASENTKSNQPEE
ELENIVGVDIPLVNEEAETTANELLDNDSEKGLTIPTLN
QSDQDNSTADASKNDESTPSPTEVHSLCTISNQETIKA
GRIPPLGERQSESMPPDRKMKNVFFPLEVKLKDNYPSIL
KADSSLSEVAFSPKEYNSPSFRPEYIEFKPKFQDCSDK
FEDLHDMTSFTHLKKVKSRSRLLGKSSDDIHNHARHS
ARPYTAPEVKNQRESYSGKFTSRMVSSGLVHINDKT
SDYEMHKMRPKKIKRGY

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 151254

Gene Symbol: ALS2CR11

Gene Alias: FLJ25351, FLJ40332